

## Q1

## Q6

Model name		SHC-U8S	SHC-U7D
Number of channels		6	1
Size		W 145mm, L 125mm, H 12mm	W 147mm, L 33mm, H 12mm
Weight		45g	22g
Color		White	White
IP grade		IP22	Not available
Measuring Signals		Sampling rate (per second)	250
		Filter range (Hz)	0.5~40
		Continuous recording (hours)	120
Power	Device	Power consumption	25mW
		Battery	Lithium-Polymer 3.8VDC, 960mAh
	Charging	Input	5.0VDC, 1A
Communication		Data transfer	Bluetooth 5.0
		Transmission distance (m)	<10
Interworking software		App	QUET ECG (Android, iOS - coming soon)
		PC	Quet PRO
		File format	RPD, CSV
Components		1 Refit Patch U8 (SHC-U8), 1 Charging cradle, 3 Clips, 3 cables	1 Refit Patch U7D (SHC-U7D), 1 Charging cradle, 2 Clips, 1 cable
Manufacturer/Country		Solmitech Co., Ltd. / Republic of Korea	Solmitech Co., Ltd. / Republic of Korea