

Device Data Sheet

Model name		SHC-U7D	
Number of channels		1	
Size		W 147mm, L 33mm, H 12mm	
Weight		22g	
Color		White	
IP grade		Not available	
Measuring Signals		Sampling rate (per second)	250
		Filter range (Hz)	0.5~40
		Continuous recording (hours)	10
Power	Device	Power consumption	37mW
		Battery	Lithium-Polymer 3.7VDC, 180mAh
	Charging	Input	5.0VDC, 1A
Communication		Data transfer	Bluetooth 2.1 / 4.0
		Transmission distance (m)	<10
Interworking software		App	QUET ECG (Android, iOS)
		PC	Quet PRO
		File format	RPD, CSV

Components	1 Refit Patch U7D (SHC-U7D), 1 Charging cradle, 2 Clips, 1 cable
Manufacturer/Country	Solmitech Co., Ltd. / Republic of Korea