

Q6 Data Sheet

Model name		SHC-U8S	
Number of channels		6	
Size		W 145mm, L 125mm, H 12mm	
Weight		45g	
Color		White	
IP grade		IP22	
Measuring Signals		Sampling rate (per second)	250
		Filter range (Hz)	0.5~40
		Continuous recording (hours)	120
Power	Device	Power consumption	25mW
		Battery	Lithium-Polymer 3.8VDC, 960mAh
	Charging	Input	5.0VDC, 1A
Communication		Data transfer	Bluetooth 5.0
		Transmission distance (m)	<10
Interworing software		App	QUET ECG (Android, iOS - coming soon)
		PC	Quet PRO
		File format	RPD, CSV
Components		1 Refit Patch U8 (SHC-U8), 1 Charging cradle, 3 Clips, 3 cables	

Manufacturer/Country	Solmitech Co., Ltd. / Republic of Korea
----------------------	---